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



Radioactive hazard

Explosive hazard

Download from  
**Dreamstime.com**  
This watermarked comp image is for previewing purposes only.

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## **Ethidium bromide safety precautions**

- Ethidium bromide is a mutagen and should be handled as a hazardous chemical, so wear a laboratory coat and protective gloves to prevent skin contact when working with liquid solutions. Change gloves frequently.
- Avoid handling the powder directly to prevent inhalation exposure, work with EtBr powder or crystals in a fume hood, or work with premixed EtBr solutions or tablets.
- Wear eye protection and ensure that there is unobstructed access to an eyewash/shower unit in the work area.
- Do not eat or drink where EtBr is handled, processed, or stored to avoid its ingestion.
- Always wash hands thoroughly after handling EtBr, even if gloves are used.
- Do not dispose gel, gloves or paper towels in normal trash receptacles. Only dispose them in appropriate biohazard boxes in your lab.
- Store EtBr away from strong oxidizing agents in a cool, dry place, and keep the container undamaged and tightly closed.

### **Emergency exposure procedures**


- **Eye care**: If EtBr comes in contact with the eyes, immediately flush them with copious amounts of cold or cool water for at least 15 minutes.
- **Skin care**: In the event of skin exposure, remove contaminated clothing and immediately wash the affected area with soap and copious amounts of cold or cool water for 15 minutes.
- **If swallowed or inhaled**: In the case of EtBr ingestion, obtain medical attention immediately. If EtBr dust is inhaled, move the victim to a source of fresh air.

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**Model: PURI STER, Serial No. PSA60-1307085, CRYSTE**

**Instructions:**

- 1) Fill water in the base of the autoclave chamber.
- 2) Place the loaded basket drum in the autoclave.
- 3) Place the items you want to sterilize them in the device.
- 4) Put on the lid, press start button.
- 5) Let the process go on and after 20 min, turn off the power when the pressure become zero lb psi.
- 6) Now open the lid and remove the sterilized items to preserve them in a safe place.
- 7) Remove all water from the base of the autoclave.

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**Model: PS.3A Advanced Technology**

**Instructions:**

- 1) Make sure that all the glass wares are dry.
- 2) Put the items you want to sterilize them after being rapped in aluminum foil/ paper.
- 3) Adjust the temperature at 180 degree Celsius.
- 4) Turn on the oven and leave the cycle for 2 h to ensure proper sterilization.
- 5) Turn off the device.

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**Model: CxL**

**Instructions:**

- 1) Place the microscope slide on the stage and fix it with the stage clips.
- 2) Rotate the focus knob to move the stage upward without letting the objective lens touch the cover slide.
- 3) Look through the eyepiece and move the focus knob until the image comes into focus.
- 4) Adjust the condenser and light intensity.
- 5) Move the slide around until the sample is in the center of the field.
- 6) Remove the slide and turn the device off after cleaning the stage and the microscope with xylene.

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**Model :** D3165 Hanigsen/ W. Germany

**Instructions**

- 1- Set the required temperature using the temperature control and leave the **Water bath** to run for 30 -60 mins.
- 2- Place a thermometer in the centre of the **Water bath**
- 3- Take the temperature reading after 30 minutes; if the temperature is not as set  $\pm 1$  °C, adjust the control button (softy button) until the required temperature is reached.
- 4- Repeat the process every 30 minutes
- 5-Use the water bath with caution to avoid the electrical hazards and burns
- 6- Clean the water bath from inside and change the water before use to prevent deposition of salts and prevent the microbial growth
- 7- Switch off the equipment at the end of the operation and when not in use.

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**Name of the instrument:** Cooling centrifuge

**Model:** PRO-RESEARCH.K2015R

**Instructions:**

- 1- Press the book sign to set a program. 
- 2- Press (**set parameters**) then select the required rotor, time and RPM.
- 3- Press tick (✓) after selection.
- 4- Press central button on right to open the lid. 
- 5- Inspect centrifuge tubes for cracks before use.
- 6- Clearly label tubes for identification.
- 7- Make a counterbalance for the centrifuge tubes you want to put in the centrifuge then close the lid pressing both sides down firmly.
- 8- Press start button. 

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**Name of the instrument:** Thermocycler

**Model:** TC96K, Serial No. GM-02-1122, Acculab

**Instructions:**

- 1) Turn on the operating button.
- 2) Load the reactions into 0.2 ml PCR tubes (use Eppendorf tubes to ensure good fit in machine).
- 3) Make sure to use small holes in the PCR machine.
- 4) You can load up to 96 tubes at a time in this device.
- 5) Close the lid well and rotate the knob to be tightly closed.
- 6) Select **File** from the list on the screen of the device.
- 7) Select your program then press **Enter**.
- 8) Select F2 (Edit) to adjust the number of segments, steps and cycles as you want according to your protocol.
- 9) Press F4 (Next) then **save**.
- 10) Press **Run** to start the process.
- 11) When the run ends, press **Stop** then **Enter**.
- 12) Remove your tubes and turn of the machine after cleaning it.

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
<b>Name of the instrument: Gel electrophoresis</b>
<b>Model: CS-300v</b>
<b>Instruction</b>
<p>1- Prepare TBE buffer and add a sufficient amount to a conical flask.</p> <p>2- Weigh agarose powder according to your protocol and add it to the buffer solution .</p> <p>3- Keep the solution-containing flask on the hot plate until the agarose powder is completely dissolved in the buffer.</p> <p>4- Allow the solution to cool for few minutes then, add ethidium bromide solution and mix well.</p> <p><b>Caution :</b> Ethidium bromide is a mutagen and should be handled as a hazardous chemical so avoid direct contact and wear protective gloves while handling.</p> <p>5- Make sure that the gel tray and comb were properly placed together. Pour the solution to the gel caster.</p> <p>6- Allow the gel to set completely at room temperature then, carefully remove the comb and place the gel in the electrophoresis tank.</p> <p>7- Add enough TBE buffer to the tank and ensure that the gel is completely immersed in the buffer.</p> <p>8- Mix the samples of DNA with the desired 6x loading dye and load the samples carefully inside the pre-casted wells in an organized manner using Automatic pipette/ add the ladder.</p> <p>9- Close the lid of the electrophoresis tank, connect the instrument to the power supply and allow it to run at 250v for 30-45 minutes.</p> <p>10- After finishing the run, switch off the power supply, remove the gel carefully to UV transilluminator and follow its instructions.</p>

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<b>Name of the instrument: Gel documentation system</b>
<b>Model : GelDoc 1.4</b>
<b>Instructions</b>
<p>1- Open up the sliding door at the bottom of the and place the prepared gel onto centre of open tray.</p> <p>2- Close the sliding door completely and turn on Gel documentation system (UV light) and the computer, select the Imaging icon then, select the “<b>Image capture</b>” button on the right side of the opened window to take photo for gel.</p> <p>3- You can zoom the camera out and manipulate the resolution of the image using the buttons on the upper side of the gel documentation instrument.</p> <p>4- If the gel is moved, close the UV transilluminator, open the door and move the gel manually and close the door after manipulation then, turn on the instrument.</p> <p>5- After Image capture, click <b>File, Save as</b>, then chose your document for save then, press <b>Save</b>.</p> <p>6- After finishing, take off gel and dispose it in the biohazard box.</p> <p>7- Do not dispose of any gels in normal trash receptacles. Only dispose of gels, gloves, paper towels in appropriate biohazard boxes in your lab or the biohazard box located in this room.</p> <p>8- Make sure you clean all equipment surfaces and table before you leave room and make sure you have shut down Gel documentation system , shut down the computer, the room lightening and door after you finish.</p> <p><b>Attention:</b></p> <ul style="list-style-type: none"><li>- Do not cut gels or use any sharp instruments on glass surface of UV transilluminator. This will damage glass surface reducing the equipment’s ability to take high quality images.</li><li>- The whole system may be contaminated with Ethidium Bromide that is a potential carcinogen so use the system with caution.</li></ul>

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<b>Name of the instrument:</b> SDS Polyacrylamide gel electrophoresis
<b>Model:</b> 150507 696
<b>Instructions:</b> <ol style="list-style-type: none"><li>1) Gather combs, glass plates, spacer and blinder clips</li><li>2) Thoroughly clean the glass plates with ethanol and assemble the gel casting mold</li><li>3) Pour acrylamide solution for a separating gel, overlay with isopropanol to prevent contact with air which inhibits polymerization.</li><li>4) Allow acrylamide to polymerize for 30 min to form a gel and then get rid of the overlaid isopropanol then wash with water to remove the traces.</li><li>5) Pour acrylamide solution for a stacking gel; insert a comb and allow the acrylamide to polymerize.</li><li>6) Prepare the samples (add sample buffer to samples and mix them well and finally heat them at 100 degree Celsius for 10 min).</li><li>7) Remove the binder clips, spacer and comb from the gel and mount the gel in the electrophoresis chamber.</li><li>8) Pour running buffer into the upper and lower chambers and remove all air bubbles.</li><li>9) Load samples and the protein marker in wells.</li><li>10) Turn on the power supply and run the gel at 250 V for 1 h until the dye in the sample buffer reaches the bottom of the gel.</li><li>11) Remove the gel assembly from the electrophoresis apparatus.</li><li>12) Remove the gel assembly from the glass plates using spatula and prepare for subsequent analysis.</li></ol>

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**SDS PAGE device**

**Name of the instrument: Western blot device**

**Model: REL-SD10**

**Instructions:**

- 1- Wear gloves and cut the nitrocellulose membrane and 4 Whatman filter papers according to the size of the gel.
- 2- Incubate the gel, membrane and the filter papers in blot buffer for two minutes.
- 3- Place two filter papers toward the anode side (down) of the semi-dry blotter, then put the membrane, the gel and the other two filter papers, respectively.
- 4- Close the lid and start the transfer at 250 mA and 130 V for 1hr.
- 5- Block the membrane by incubation in 5% skimmed milk in TBS-T buffer at room temperature for 1 hr.
- 6- For His-tagged proteins: incubate the blocked membrane with anti-His-HRP Ab in TBS-T for 1 hr with gentle shaking at room temperature.
- 7- Wash the membrane with TBS-T buffer for 3\*5 min.
- 8- Expose the membrane to TMB substrate solution to detect the protein bands.

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
**Name of the instrument: BioTek ELISA reader**

**Model: VT 05404-0998**

**Instructions:**

- 1- Switch on the BioTek ELISA reader and the computer device.
- 2- Open Gen5 software, press **read now**→ **New**.
- 3- If you are using lid press  $\checkmark$  beside (**use lid**).
- 4- Press (**Read**) → select the wavelength you want from the inserted wavelength options (**405, 450, 490, 630**).
- 5- If you don't want to read the whole plate press (**full plate**) box and select the cells you want to read.
- 6- Press **Ok** to start reading.
- 7- Save your data from **file** →**save as** →**excel sheet**.

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© Jozsef Szasz-fabian | Dreamstime.com